

EVERPURE® MICROGUARD HF POST-FILTER

PART NO. PWMG13



Microguard HF Post-Filter Part no. PWMG13

APPLICATIONS

Small capacity water coolers and drinking fountains.

Use as a post-filter at or near dispense tap.

For use with still drinking water up to 40°C.

OPERATION TIPS

Change cartridges on a regular 6 month preventative maintenance program or when flow rate drops significantly.

Always use prefiltration to maximise membrane life.

On installation, flush filter at full flow for 2 minutes.

FEATURES • BENEFITS


3-inch capillary membrane inline filter with micro-filtration (MF) membranes.


MF reduces 99.9999% of bacteria such as Pseudomonas, cysts such as Cryptosporidium and Giardia as well as asbestos fibres.

SIZING

Flow rate: 3,2 L/min @ 2 bar

Service life: Depends on local water conditions and level of pre-filtration. European guidelines recommend exchange of cartridge after 6 months.





SPECIFICATIONS

Membrane tupe:

Micro-filtration membranes; asymmetric, hydrophilic capillaries

Membrane material:

PES/PVP

Max. pore size:

0,15µm absolute

Membrane surface:

Approx. 280cm²

Retention:

Proven>99,9999% (>log 6) of bacteria and cysts according to; ASTM F838-05

Material:

Housing and connectors: PS Swivel and

filter tube:PVC Potting: 2K epoxy resin

Max. pressure:

6 bar (88 psi) 40°C (140°F)

Max. temp:

Cartridge size:

33 x 210 mm (OD x Length)


Cartridge:

Cartridge tested for integrity according

to ASTM D6908-03.

Weight:

28g





Liberty Ersatzteile & Zubehör

Liberty Vertriebs GmbH Bonner Straße 103 40589 Düsseldorf

Tel.: 0211-753520 Fax.: 0211-753575 Internet: www.libertyshop.de Email: verkauf@libertyshop.de